

rS Red

Seat Colors Charcoal/Gray Vinyl

Body Colors
Blue Red



The heavy-duty version of the original Compass offers users a robust 450 lb. weight capacity with all of the same great features including the patented springless articulating chassis and true center-wheel drive design for unmatched maneuverability in tight spaces. The Compass HD comes in a variety of seat sizes up to 24" wide.



YOUR LIFE IN MOTION



Golden Compass[™] HD Specifications

dolucii dollipass	TID Opcome	autons
Model	GP620CC	GP620SS
Model Name	Compass HD	Compass HD
Medicare Code	K0825	K0824
Weight Capacity	450 lbs	450 lbs
Drive Wheels	Center	Center
Maximum Speed	4 mph	4 mph
Operating Range ¹	12 Miles	12 Miles
Ground Clearance	3.5"	3.5"
Turning Radius	20"	20"
Type Batteries	2-NF22	2-NF22
Colors	Red. Blue	Red, Blue
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes
Charger	Off Board	Off Board
•	Dynamic 75 Amp	Dynamic 75 Amp
Controller	Remote	Remote
Dimensions:		
Length	40"	40"
Width	25"	25"
Height (ground to top of back of seat)	38.25" - 42.75"	38.25" - 42.75"
Height (ground to top of back of headrest)	47.25" - 53.75"	47.25" - 53.75"
Ground to Top of Seat	21.25" - 24.25"	19.75" - 22.25"
Number of Seat Height Positions	5	5
Size of Increments	0.5"	0.5"
Weight of Unit: (assembled)	226 lbs	223 lbs
Base	101 lbs	101 lbs
Seat w/Arms	51 lbs	48 lbs
Batteries	37 lbs ea	37 lbs ea
Tires:	Flat Free Black	Flat Free Black
Center	10"	10"
Casters	6"	6"
Standard Seat:	Captain's Highback	Captain's Van/Pan
Back Height - no headrest	21"	23"
Back Height with headrest	28.5" - 32.5"	30.5" - 34.5"
Width x Depth	20" X 20", 22" X 20", 24" X 20"	20" X 20"
Color	Coffee/Sand, Charcoal/Gray	Coffee/Sand, Charcoal/Gray
Optional Seats:	Captain's Highback	Captain's Van/Pan
Back Height - no headrest	21"	23"
Back Height with headrest	28.5" - 32.5"	30.5" - 34.5"
Width	18", 20", 22" or 24"	18", 20", 22" or 24"
Depth	18" or 20"	16" or 18"
Color	Coffee/Sand, Charcoal/Gray	Coffee/Sand, Charcoal/Gray
Optional Seats 2:	N/A	Pro-Flex Seat
Width	N/A	18", 20" or 22"
Depth	N/A	18" or 20"
Color	N/A	Black
Footrest:	Yes	Yes
Height adjustable	3" - 7"	3" - 7"
Angle adjustable	0 - 35 degrees	0 - 35 degrees
Warranty:		
Frame	5 Years	5 Years
Drive Train	13 Months	13 Months
Electronics ²	13 Months	13 Months

¹ Battery range will vary due to rider weight, drive surface and terrain.

^{2.} Electronics warranty excludes batteries.



Power Chair Accessories and Options

- Captain's seat (GP620CC)
- Pro Flex seat (GP620SS)
- Van/pan seat (GP620SS)
- Swing-away legrests
- Elevating legrests
- Stump support
- Footplate extender
- Cane/crutch holder kit
- Walker holder kit
- Oxygen tank holder kit
- Pack N Go

2 exciting colors





Red

Blue

