Profile Mid-Wheel Drive Lift AL580 / AL580XL





Specifications	AL580	AL580XL Extra Large
Accommodates	Mid-wheel power chairs with 10" drive wheels	Mid-wheel power chairs with 14" drive wheels
Platform	16.5" x 27"	17.5" x 27"
Lifting Capacity	350 lbs. (160 kg)	350 lbs. (160 kg)
Installed Weight	63 lbs. (31 kg)	63 lbs. (31 kg)
Securement	Padded hold-down arm & wheel cradle	Padded hold-down arm & wheel cradle
Hitch	Class II or greater	Class II or greater
Included Accessories	AL123 Class III hitch adapter / AL107 hitch post	AL123 Class III hitch adapter / AL107 hitch post
Options	Battery pack, Swing-Away, Power chair covers	Battery pack, Swing-Away, Power chair covers
Warranty	3-year transferable	3-year transferable







