## **Options**

- Distance straps A variety of Human Care's distance straps in different lengths can be used to achieve a suitable lifting height for both the patient and the caregiver. Must allow minimal distance of 72.8" (1850 mm) between the top position of the bar and the floor
- Four-point suspension Attaches to lift's single leg hook for versatile and comfortable patient suspension. Can be used with stretcher accessory for horizontal lifting
- Human Care Data and Service Application Collects data and information about the lift's usage to ensure proper use, minimize risk of back injury and increase return on equipment investment. Can access information from PC or tablet device
- **Lift trolley** Allows for easy storage and transportation of Altair, its slings and other corresponding accessories
- Niklas freestand Altair attaches onto this aluminum gantry to provide a safe and easy moving method for caregivers. This freestanding system requires no assembly or modifications to walls or ceilings
  - Weight capacity: 485 lbs (220 kg)
  - Unit dimensions:
    - Minimum: 62" x 63.5" x 45.1" (1579 mm x 1615 mm x 1145 mm)
    - Maximum: 94" x 103.2" x 45.1" (2388 mm x 2620 mm x 1145 mm)
- Reacher arm Lightweight metal arm aids user in hanging up and taking down Altair. Arm length: 23" (584 mm)
- Scale Medical grade scale can weigh patients up to 705.5 lbs (320 kg)
- Slings A variety of Human Care's slings are compatible with its overhead lift motors to achieve safe, versatile and comfortable patient suspension and transfer. With over 22 different sling models, sized from junior small to bariatric, Human Care slings can be used for:
  - General lifting situations
  - Hygiene lifting in connection to toiling or bathing
  - Repositioning
  - o Gait training and rehabilitation
  - Bariatric lifting
  - Specialized lifting situations

Additionally, Human Care slings are available in three different materials, including:

- o Polyester
- Polyester net
- Trevira

- **Stretcher** Light, portable stretcher enables easy horizontal lifting. Must be used with four-point suspension bar.
  - o Weight capacity: 441 lbs (200 kg)
  - Unit dimensions: 80.3" x 29.1" x 1.6" (2040 mm x 740 mm x 40 mm)
- **Weight capacity** Three different weight capacities for optimum customization:
  - o 331 lbs (150 kg)
  - o 485 lbs (220 kg)
  - o 661 lbs (300 kg)