



True, stable platform, center-wheel drive power chair with an articulating chassis for the ultimate in maneuverability and stability.



YOUR LIFE IN MOTION



Golden Compass™ Specifications

Model:	GP-615 CC	GP-615SS	GP-615SS
Weight:			(specs with tilt system)
Base	90 lbs	90 lbs	90 lbs
Seat	35 lbs	33 lbs	115 lbs w/armrests & electronics
Armrests	14 lbs	14 lbs	
Total Weight	217 lbs (w/ 22NF batt.)	215 lbs (w/ 22NF batt.)	283 lbs (w/ 22NF batt.)
Overall Dimensions:			
Length	39.8"	39.8"	42.7"
Width	25"	25"	26.2"
Seat Height from ground	20.5" (w/ 22NF batt.) adjustable to 23.5"	19" (w/ 22NF batt.) adjustable to 22"	20" (w/ 22NF batt.) adjustable to 23"
Top Speed:	6.5 mph	6.5 mph	6.5 mph
Operating Range:	12 miles*	12 miles*	12 miles*
Load Capacity:	300 lbs	300 lbs	300 lbs
Turning Radius:	20"	20"	21.95"
Wheel Sizes:			
Drive Wheels	10" x 3.5"	10" x 3.5"	10" x 3.5"
Casters	6" x 2"	6" x 2"	6" x 2"
Seat Dimensions:			
Standard Captain's:	18" x 18"	18" x 18"	18" x 18"
Max. dimension	24" x 20"	24" x 20"	22" x 20"
Min. dimension	16"x 16"	16"x 16"	14"x 14"
Footrest:			
Height adjustable	Yes	Yes	Yes
Angle adjustable	Yes	Yes	Yes
Length adjustable	Yes	Yes	Yes



Controller: Shark 75 amp remote controller, fully programmable GP-615 Standard: **Charger:** Off Board Charger Colors: Candy Apple Red, Blue, Onyx, Hunter Green **Seat Colors:** Standard: Coffee/Sand, Charcoal/Gray Captain's Optional: Black Vinyl, Stadium Style

Battery Type (recommended): 2 (two) 12 Volt 22NF

Power Chair Accessories and Options

- Stadium seat
- High back seat
- Pro Flex seat Van/pan seat

- Elevating legrests
- Stump support
- Footplate extender
- Positioning belt
- Swing-away legrests
 Cane/crutch holder kit
 Pack N Go
- Walker holder kit
- Oxygen tank holder kit



Compass TRO shown with Motion Concepts' Power Tilt

4 exciting colors



Candy Apple Red



Blue



Hunter Green





^{*} Battery range will vary due to rider weight, drive surface, drive slope Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.