


- Extended hook to clear rear bumper
- Red straps included for tailgate usage
- Separates into two pieces for easy carrying
- Rotation of pin either locks or allows ramp to be separated into two pieces
- Folds, carries like a suitcase
- Patented hinge eliminates pinch point
- Closure strap locks ramp panels together
- Lightweight, easy-to-handle & set up
- Durable welded construction
- Full platform provides excellent stability
- Anti-slip, high traction surface
- Accommodates wheelchairs and scooters with various wheel configurations
- Safety DVD and steel security pins included



Safety guide indicates if ramp is on a safe slope!

Item Number	Length	Width	Folded Size	Weight	Capacity*
UTW630	6'	30"	38 $\frac{1}{4}$ " x 16" x 10 $\frac{3}{4}$ "	36 lbs.	800 lbs.
UTW730	7'	30"	44 $\frac{1}{4}$ " x 16" x 10 $\frac{3}{4}$ "	41 lbs.	800 lbs.
UTW830	8'	30"	50 $\frac{3}{4}$ " x 16" x 11"	52 lbs.	800 lbs.
UTW1030	10'	30"	62 $\frac{3}{4}$ " x 16" x 11"	63 lbs.	700 lbs.

Not sure which ramp length you need? Refer to pages 18 - 19 to help you.

LIMITED LIFETIME WARRANTY

* Maximum ramp capacity: 400 lbs. (one axle) – 800 lbs. (two axles). Except on UTW1030: 350 lbs. (one axle) – 700 lbs (two axles).

